

Trade name: tri-Cobalt tetra-oxide**Product no.:** 638-60239-41**Current version :** 4.0.0, issued: 04.01.2022**Replaced version:** 3.2.0, issued: 02.11.2020**Region:** GB**SECTION 1: Identification of the substance/mixture and of the company/undertaking****1.1 Product identifier****Trade name****tri-Cobalt tetra-oxide**

Substance name tricobalt-tetraoxide

Identification numbers

CAS no. 1308-06-1

EC no. 215-157-2

1.2 Relevant identified uses of the substance or mixture and uses advised against**Relevant identified uses of the substance or mixture**

General chemical processes

Uses advised against

No data available.

1.3 Details of the supplier of the safety data sheet**Address**

Shimadzu Europa GmbH

Albert-Hahn-Str. 6-10

47269 Duisburg

Telephone no. +49 (0)203 - 7687-0

Fax no. +49 (0)203 - 711045

Information provided by / telephone

Tel.: +49 (0)203 - 7687-0

Advice on Safety Data Sheet

sdb_info@umco.de

1.4 Emergency telephone number

For medical advice (in German and English):

+49 (0)551 192 40 (Giftinformationszentrum Nord)

SECTION 2: Hazards identification**2.1 Classification of the substance or mixture****Classification in accordance with Regulation (EC) No 1272/2008 (CLP)**

Carc. 1A; H350i

Carc. 2; H351

Skin Sens. 1; H317

Classification information

This product is assessed and classified using the methods and criteria below referred to in Article 9 of Regulation (EC) n° 1272/2008:

Physical hazards: determined through assessment data based on the methods or standards referred to in part 2 of Annex I to CLP

Health hazards and environmental hazards: determined through toxicological and ecotoxicological assessment data based on the methods or standards referred to in Part 3 and 4 of Annex I to CLP.

2.2 Label elements**Labelling according to Regulation (EC) No 1272/2008 (CLP Regulation)****Product identifier**

1308-06-1 (tricobalt-tetraoxide)

Hazard pictograms

GHS07



GHS08

Signal word

Danger

Hazard statement(s)

H317

May cause an allergic skin reaction.

H350i

May cause cancer by inhalation.

H351

Suspected of causing cancer.

Precautionary statement(s)

P201

Obtain special instructions before use.

Trade name: tri-Cobalt tetra-oxide**Product no.:** 638-60239-41**Current version :** 4.0.0, issued: 04.01.2022**Replaced version:** 3.2.0, issued: 02.11.2020**Region:** GB

P261 Avoid breathing dust.
 P280 Wear protective gloves/protective clothing.
 P308+P313 IF exposed or concerned: Get medical advice/attention.
 P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

Supplemental label elements
 "Restricted to professional users"

2.3 Other hazards

No data available.

SECTION 3: Composition/information on ingredients**3.1 Substances****Chemical characterization**

Substance name tricobalt-tetraoxide
 Formula Co₃O₄
 Molecular weight 240.8

Identification numbers

CAS no. 1308-06-1
 EC no. 215-157-2

Acute toxicity estimate (ATE) values

oral	dermal	inhalative
		4,84 mg/l

3.2 Mixtures

Not applicable. The product is not a mixture.

SECTION 4: First aid measures**4.1 Description of first aid measures****General information**

In case of persisting adverse effects, consult a physician. Remove contaminated clothing and shoes immediately, and launder thoroughly before reusing.

After inhalation

Remove affected person from the immediate area. Ensure supply of fresh air. Keep breathing passages free. Irregular breathing/no breathing: artificial respiration. Call a doctor immediately.

After skin contact

Wash off immediately with soap and water. Seek medical attention.

After eye contact

Separate eyelids, wash the eyes thoroughly with water (15 min.). Seek medical assistance.

After ingestion

Call a doctor immediately. Never give anything by mouth to an unconscious person. If individual is drowsy or unconscious, place in recovery position (on left side, with head down). Rinse the mouth thoroughly with water.

4.2 Most important symptoms and effects, both acute and delayed

No data available.

4.3 Indication of any immediate medical attention and special treatment needed

No data available.

SECTION 5: Firefighting measures**5.1 Extinguishing media****Suitable extinguishing media**

Water; Carbon dioxide; Extinguishing powder

Unsuitable extinguishing media

No data available.

5.2 Special hazards arising from the substance or mixture

In the event of fire, the following can be released: Cobalt dust

5.3 Advice for firefighters

Use self-contained breathing apparatus. Wear protective clothing.

SECTION 6: Accidental release measures**6.1 Personal precautions, protective equipment and emergency procedures**

Trade name: tri-Cobalt tetra-oxide**Product no.:** 638-60239-41**Current version :** 4.0.0, issued: 04.01.2022**Replaced version:** 3.2.0, issued: 02.11.2020**Region:** GB**For non-emergency personnel**

Refer to protective measures listed in sections 7 and 8. Avoid dust formation. Keep away from ignition sources. Ensure adequate ventilation. Avoid contact with skin, eyes and clothing.

For emergency responders

No data available. Personal protective equipment (PPE) - see Section 8.

6.2 Environmental precautions

Do not discharge into the drains/surface waters/groundwater.

6.3 Methods and material for containment and cleaning up

Collect mechanically. Send in suitable containers for recovery or disposal.

6.4 Reference to other sections

No data available.

SECTION 7: Handling and storage**7.1 Precautions for safe handling****Advice on safe handling**

Provide good ventilation at the work area (local exhaust ventilation, if necessary). If workplace exposure limits are exceeded, respiratory protection approved for this particular job must be worn. Avoid the formation and deposition of dust.

General protective and hygiene measures

Avoid contact with eyes and skin. Do not inhale dust. At work do not eat, drink, smoke or take drugs. Keep away from foodstuffs and beverages. Remove soiled or soaked clothing immediately. Clean skin thoroughly after work; apply skin cream. Store work clothing separately. Provide eye wash fountain in work area. Have emergency shower available.

Advice on protection against fire and explosion

Dust can form an explosive mixture with air. Take precautionary measures against static charges. Keep away from sources of heat and ignition.

7.2 Conditions for safe storage, including any incompatibilities**Technical measures and storage conditions**

Keep container tightly closed and dry in a cool, well-ventilated place.

Requirements for storage rooms and vessels

Always keep in containers of same material as the original. Store product in closed containers.

Incompatible products

None known

7.3 Specific end use(s)

No data available.

SECTION 8: Exposure controls/personal protection**8.1 Control parameters****Occupational exposure limit values**

No	Substance name	CAS no.	EC no.
1	tricobalt-tetraoxide	1308-06-1	215-157-2
List of approved workplace exposure limits (WELs) / EH40			
Cobalt & cobalt compounds (as Co)			
	WEL long-term (8-hr TWA reference period)	0.1	mg/m ³
	Comments	Carc (cobalt dichloride and sulphate), Sen	

8.2 Exposure controls**Appropriate engineering controls**

No data available.

Personal protective equipment**Respiratory protection**

If workplace exposure limits are exceeded, a respiration protection approved for this particular job must be worn. In case of dust formation, take appropriate measures for breathing protection in the event that workplace threshold values are not specified.

Eye / face protection

Tightly fitting safety glasses (EN 166).

Hand protection

Sufficient protection is given wearing suitable protective gloves checked according to i.e. EN 374, in the event of risk of skin contact with the product. Before use, the protective gloves should be tested in any case for its specific work-station suitability (i.e. mechanical resistance, product compatibility and antistatic properties). Adhere to the manufacturer's instructions and information relating to the use, storage, care and replacement of protective gloves. Protective gloves shall be replaced immediately when physically damaged or worn.

Trade name: tri-Cobalt tetra-oxide**Product no.:** 638-60239-41**Current version :** 4.0.0, issued: 04.01.2022**Replaced version:** 3.2.0, issued: 02.11.2020**Region:** GB

Design operations thus to avoid permanent use of protective gloves.

Other

Normal chemical work clothing.

Environmental exposure controls

No data available.

SECTION 9: Physical and chemical properties**9.1 Information on basic physical and chemical properties**

State of aggregation	
solid	
Form/Colour	
granules	
black	
Odour	
odourless	
pH value	
No data available	
Boiling point / boiling range	
No data available	
Melting point/freezing point	
Value	895 °C
Comments	The substance decompose by heating.
Decomposition temperature	
Value	> 900 °C
Flash point	
Not applicable	
Ignition temperature	
No data available	
Flammability	
No data available	
Lower explosion limit	
No data available	
Upper explosion limit	
No data available	
Vapour pressure	
No data available	
Relative vapour density	
No data available	
Relative density	
No data available	
Density	
Value	6.07 g/cm ³
Solubility in water	
Comments	insoluble
Solubility	
No data available	
Partition coefficient n-octanol/water (log value)	
No data available	
Viscosity	
No data available	
Particle characteristics	
No data available	

9.2 Other information

Other information
No data available.

Trade name: tri-Cobalt tetra-oxide

Product no.: 638-60239-41

Current version : 4.0.0, issued: 04.01.2022

Replaced version: 3.2.0, issued: 02.11.2020

Region: GB

SECTION 10: Stability and reactivity**10.1 Reactivity**

Dangerous reactions are not expected if the product is handled according to its intended use.

10.2 Chemical stability

Stable under recommended storage and handling conditions (See section 7).

10.3 Possibility of hazardous reactions

No data available.

10.4 Conditions to avoid

Protect from humid air and water.

10.5 Incompatible materials

Reducing agents; Acids

10.6 Hazardous decomposition products

Cobalt oxide

SECTION 11: Toxicological information**11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008**

Acute oral toxicity			
No	Substance name	CAS no.	EC no.
1	tricobalt-tetraoxide	1308-06-1	215-157-2
LD50	>	5000	mg/kg bodyweight
Species	rat		
Acute dermal toxicity			
No	Substance name	CAS no.	EC no.
1	tricobalt-tetraoxide	1308-06-1	215-157-2
LD50	>	2000	mg/kg bodyweight
Species	rat		
Acute inhalational toxicity			
No	Product Name	CAS no.	EC no.
LC50		4.83	mg/l
Duration of exposure		4	h
Species	rat		
No	Substance name	CAS no.	EC no.
1	tricobalt-tetraoxide	1308-06-1	215-157-2
LC50		4.83	mg/l
Duration of exposure		4	h
State of aggregation	Dust		
Species	rat		
Skin corrosion/irritation			
No data available			
Serious eye damage/irritation			
No data available			
Respiratory or skin sensitisation			
No data available			
Germ cell mutagenicity			
No data available			
Reproduction toxicity			
No data available			
Carcinogenicity			
No data available			
STOT - single exposure			
No data available			
STOT - repeated exposure			
No data available			
Aspiration hazard			
No data available			
Delayed and immediate effects as well as chronic effects from short and long-term exposure			

Trade name: tri-Cobalt tetra-oxide**Product no.:** 638-60239-41**Current version :** 4.0.0, issued: 04.01.2022**Replaced version:** 3.2.0, issued: 02.11.2020**Region:** GB

Cobalt and cobalt compounds may lead in case of exposition to higher concentrations to symptoms of poisoning and may cause damage to kidneys, the nervous and cardiovascular system. Allergical reactions like skin eczema possible. Sensitisation possible through inhalation. Inhalation of dusts may irritate the respiratory tract. Frequent persistent contact with the skin can cause skin irritation.

11.2 Information on other hazards**Endocrine disrupting properties**

No data available.

Other information

No data available.

SECTION 12: Ecological information**12.1 Toxicity****Toxicity to fish (acute)**

No data available

Toxicity to fish (chronic)

No data available

Toxicity to Daphnia (acute)

No data available

Toxicity to Daphnia (chronic)

No data available

Toxicity to algae (acute)

No data available

Toxicity to algae (chronic)

No data available

Bacteria toxicity

No data available

12.2 Persistence and degradability

No data available.

12.3 Bioaccumulative potential

No data available.

12.4 Mobility in soil

No data available.

12.5 Results of PBT and vPvB assessment

No data available.

12.6 Endocrine disrupting properties

No data available.

12.7 Other adverse effects

No data available.

12.8 Other information**Other information**

Do not discharge into drains or waters and do not dispose of in public landfills.
 Ecological data are not available.

SECTION 13: Disposal considerations**13.1 Waste treatment methods****Product**

Allocation of a waste code number, according to the European Waste Catalogue, should be carried out in agreement with the regional waste disposal company.

Packaging

Contaminated packaging should be emptied as far as possible and after appropriate cleansing may be taken for reuse. Packaging that cannot be cleaned should be disposed of in agreement with the regional waste disposal company.

SECTION 14: Transport information**14.1 Transport ADR/RID/ADN**

The product is not subject to ADR/RID/ADN regulations.

14.2 Transport IMDG

Trade name: tri-Cobalt tetra-oxide**Product no.:** 638-60239-41**Current version :** 4.0.0, issued: 04.01.2022**Replaced version:** 3.2.0, issued: 02.11.2020**Region:** GB

The product is not subject to IMDG regulations.

14.3 Transport ICAO-TI / IATA

The product is not subject to ICAO-TI / IATA regulations.

14.4 Other information

No data available.

14.5 Environmental hazards

Information on environmental hazards, if relevant, please see 14.1 - 14.3.

14.6 Special precautions for user

No data available.

14.7 Maritime transport in bulk according to IMO instruments

Not relevant

SECTION 15: Regulatory information**15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture****EU regulations****Regulation (EC) No 1907/2006 (REACH) Annex XIV (List of substances subject to authorisation)**

In accordance with the REACH regulation (EC) 1907/2006, the product does not contain any substances that are considered as subject to listing in annex XIV, inventory of substances requiring authorisation.

REACH candidate list of substances of very high concern (SVHC) for authorisation

In accordance with article 57 and article 59 of the Reach regulation (EC) 1907/2006, this substance is not considered as subject to listing in annex XIV, inventory of substances requiring authorisation ("Authorization list").

Regulation (EC) No 1907/2006 (REACH) Annex XVII: RESTRICTIONS ON THE MANUFACTURE, PLACING ON THE MARKET AND USE OF CERTAIN DANGEROUS SUBSTANCES, MIXTURES AND ARTICLES

The substance is not subject to the provisions of annex XVII (restriction entries) of the Reach regulation (EC) 1907/2006.

Directive 2012/18/EU on the control of major-accident hazards involving dangerous substances

This substance is not subject to Part 1 or 2 of Annex I

Other regulations

Adhere to the national sanitary and occupational safety regulations when using this product.

15.2 Chemical safety assessment

No data available.

SECTION 16: Other information**Sources of key data used to compile the data sheet:**

Regulation (EC) No 1907/2006 (REACH), 1272/2008 (CLP) as amended in each case.

Directives 2000/39/EC, 2006/15/EC, 2009/161/EU, (EU) 2017/164.

National Threshold Limit Values of the corresponding countries as amended in each case.

Transport regulations according to ADR, RID, IMDG, IATA as amended in each case.

The data sources used to determine physical, toxic and ecotoxic data, are indicated directly in the corresponding section.

Creation of the safety data sheet

UMCO GmbH

Georg-Wilhelm-Str. 187, D-21107 Hamburg

Tel.: +49 40 / 555 546 300 Fax: +49 40 / 555 546 357 e-mail: umco@umco.de

This information is based on our present knowledge and experience.

The safety data sheet describes products with a view to safety requirements.

It does not however, constitute a guarantee for any specific product properties and shall not establish a legally valid contractual relationship.

Document protected by copyright. Alterations or reproductions require the express written permission of UMCO GmbH.

Prod-ID 642050